



SOPHISTICATION AND SIMPLICITY IN ONE PROGRAM

Liaison™ reads the sounds of your active life and intuitively adapts to keep you connected. As you move between different listening situations, Liaison analyzes the environment, classifies the listening situation and automatically enhances your listening experience by adapting these sophisticated features within a single program:

- ▶ Adaptive directional microphone system
- ▶ Intelligent noise reduction
- ▶ Realtime feedback canceller
- ▶ Wind noise manager

PRACTICAL SOLUTIONS TO EVERYDAY PROBLEMS

For 40 years, Unitron Hearing has been committed to making life better for people with hearing loss. To us, this means a commitment to not only develop high-quality hearing solutions, but to incorporate special features that really solve the everyday problems and concerns you have with hearing loss and hearing aids. To you, this means a company that's listening – and responding – to what matters in your life.

www.unitronhearing.com

liaison[™]
digital hearing system



ACTIVE.
INTELLIGENT.
AUTOMATIC.



THE AUTOMATIC CHOICE FOR PEOPLE ON THE GO

Fully automatic: move effortlessly between quiet and noisy environments without program changes.

Focused on speech: automatically suppresses unwanted sounds – even if they're moving – to keep you connected with the voices you want to hear.

Noise-wise: you identify what constitutes unwanted background noise – just once. From then on, Liaison remembers.

Unsurpassed comfort: no more whistling or uncomfortable wind noise.

Phone friendly: use the telephone or cell phone naturally, without switches or feedback.

Ready to go: instantly locate your ideal volume.

ALWAYS ON THE MOVE

Adaptive Directional Microphones

Liaison's adaptive directional microphone system suppresses unwanted noise from virtually any direction – even if it's moving. As your listening situation changes, Liaison's microphones automatically adapt to keep speech the priority. In quiet settings, Liaison amplifies sounds all around you. As your environment becomes noisier, Liaison intuitively tracks and suppresses unwanted sounds to maximize speech understanding.

HEAR WHAT YOU WANT

Intelligent Noise Reduction

Noise reduction systems improve comfort by tuning out unwanted background sounds. Because one person's noise might be someone else's natural connection to their surroundings, Liaison lets you decide. During fitting, your hearing healthcare professional will establish what levels of background noises you want to hear – and don't want to hear – and will program your noise reduction to engage accordingly.

ATTACKS WHISTLING

Realtime Feedback Canceller

For those who wear hearing instruments, everyday realities from phones to hats to hugs can cause uncomfortable whistling, known as feedback. Liaison's realtime feedback canceller activates instantly to prevent annoying whistling before it can be perceived, without affecting speech clarity.


ENJOY OUTDOOR PURSUITS

Wind Noise Manager

Liaison's wind noise manager intuitively adapts based on wind conditions. Wind noise is suppressed to preserve the quality and clarity of speech and prevent discomfort.

A full range of styles in a variety of signature colors

Liaison comes in a wide range of discreet styles, sizes and colors: talk to your hearing healthcare professional to determine the right choices for your needs.



Hearing healthcare professionals strongly recommend the use of hearing instruments in both ears. Your hearing healthcare professional can provide complete details.

EVERYDAY CONVENIENCES

Liaison offers the ideal blend of sophistication and practical features. Here are just two of Liaison's everyday conveniences.

ClearCall: ClearCall does away with the challenges of telephone use, allowing you to use the phone and cell phone without feedback or program changes.

Ideal volume indicator: No more fiddling with the volume control – when you reach your perfect setting on the volume control, Liaison emits a convenient little beep to let you know.*

** Styles equipped with manual volume control.*



LIAISON™
CONNECTING PEOPLE
WITH PEOPLE.