



THREE LISTENING PROGRAMS

Unison 6 and Unison 3 feature up to three listening programs for your different listening environments. As you move between environments – for example, reading at home to lunch with friends – simply change programs for the best listening experience.

JUST WHAT YOU WANT TO HEAR

Unison helps you enjoy conversations, even in noise, thanks to a dual directional microphone system in Unison 6 and Unison 3 that puts priority on sounds directly in front of you and suppresses those from other directions. Unison also provides pure sound quality by automatically suppressing very soft sounds that you don't want to hear, such as the hum of an air conditioner, and enhancing moderately soft sounds that you do, such as speech.

ABOUT UNITRON HEARING

For 40 years, Unitron Hearing has been committed to making life better for people with hearing loss. To us, this means a commitment to not only develop high-quality hearing solutions, but to incorporate special features that really solve the everyday problems and concerns you have with hearing loss and hearing instruments. To you, this means a company that's listening – and responding – to what matters in your life.

www.unitronhearing.com



THE DIGITAL CHOICE

IT'S ALL ABOUT CHOICE

The best choice of hearing instrument is different for everyone: it comes down to your individual hearing loss, lifestyle, preferences and budget. The Unison series meets your unique needs with three high-quality digital hearing instrument models.

Unison 6

An ideal choice if you have an active lifestyle and interact in a wide variety of environments.

Unison 3

An ideal choice if you're exposed to fewer listening environments.

Unison Essential

An ideal choice if you're concerned about budget and/or are making the transition from an analog to a digital hearing instrument.



unison™

COMFORTABLE LISTENING

Unison 6 controls whistling, or feedback, that can happen with any hearing instrument and is caused by fit and everyday occurrences, such as when you're chewing or pulling on a hat.

UNISON STYLE

Unison is available in the full range of shell styles, sizes and colors. Your hearing healthcare professional will help you choose the right model and style for your needs.

EVERYDAY CONVENIENCES

Ideal Volume Indicator – emits a single, discreet beep at your correct volume setting.

Easy t-coil – removes the need for a manual adjustment for phone use.

Battery Wise – maximizes battery life and provides advanced warning when batteries become low.

